

EXAMPLE OF LABELS PHYSICA ZUK

Label of the first code of each kind of device for Physica ZUK is following reported. For each label, all the information which changes in the other labels of the same device is showed.

The labels are reported as follows:

Section 1	Example of labels for all poly tibial plate (UHMWPE + AISI316L)
Section 2	Example of labels for femoral component (CoCrMo + PMMA)
Section 3	Example of labels for tibial baseplate (Ti6Al4V + PMMA)
Section 4	Example of labels for articular surface (UHMWPE)

STERILIZATION – DATA ON THE LABELS

The information about the sterilization method for the device is reported on the labels. The external cardboard box (first label in each of the following pages) indicates only the possible sterilization methods:



- **STERILE R:** sterilization by **irradiation** (Beta rays or Gamma rays according to the applicable sterilization method).

Always refer to device labels for information about the sterilization method.

For the information on applicable sterilization method for each device refer to Design Dossier.

SECTION 1 – EXAMPLE OF LABELS FOR ALL POLY TIBIAL PLATE (UHMWPE + AISI316L)

EDI

EDI 00584203110 REF

REF 5842-31-08

PHYSICA ZUK Unicompartmental Knee System

All-Poly Tibial Component

SIZE #1 h 8 LEFT Medial/RIGHT Lateral

2024-11-30 **LOT** 201900001 **STERILE R** 1900002



Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

For LimaCorporate Use Only 01



(20)01(240)96170501080(10)1900001(91)1900002(92)7008



(01)08033390125512(17)241130

Single use Do not use if package is damaged Keep away from sunlight temperature mit

Use By REF Reference Number Keep dry

LOT Lot Number IFU Consult instructions for use Caution **STERILE EO** Sterilized using Ethylene Oxide

Manufacturer **STERILE R** Sterilized using irradiation

REF 5842-31-08 All-Poly Tibial Component

PHYSICA ZUK UHMWPE+AISI 316L

Unicompartmental Knee System

#1 h 8

LEFT Medial/RIGHT Lateral

Femoral Size A,B,C,D,E,F,G



EN All-Poly Tibial Comp **DE** All-Poly Tibiakomp.

Cemented Use Only **IT** Comp. Tib. All-Poly **FR** Comp.Tibial All-Poly

LOT 201900001 **ES** Comp.Tibial All-Poly **PT** Comp.Tibial All-Poly

2024-11-30 **qty:** 1 **STERILE R** Limacorporate S.p.A.

EDI	SIZE	REF
00584203110	#1 H. 10mm	5842-31-10
00584203112	#1 H. 12mm	5842-31-12
00584203114	#1 H. 14mm	5842-31-14
00584204108	#2 H. 8mm	5842-41-08
00584204110	#2 H. 10mm	5842-41-10
00584204112	#2 H. 12mm	5842-41-12
00584204114	#2 H. 14mm	5842-41-14
00584205108	#3 H. 8mm	5842-51-08
00584205110	#3 H. 10mm	5842-51-10
00584205112	#3 H. 12mm	5842-51-12
00584205114	#3 H. 14mm	5842-51-14
00584206108	#4 H. 8mm	5842-61-08
00584206110	#4 H. 10mm	5842-61-10
00584206112	#4 H. 12mm	5842-61-12
00584206114	#4 H. 14mm	5842-61-14
00584207108	#5 H. 8mm	5842-71-08
00584207110	#5 H. 10mm	5842-71-10
00584207112	#5 H. 12mm	5842-71-12
00584207114	#5 H. 14mm	5842-71-14
00584208108	#6 H. 8mm	5842-81-08
00584208110	#6 H. 10mm	5842-81-10
00584208112	#6 H. 12mm	5842-81-12
00584208114	#6 H. 14mm	5842-81-14

EDI

EDI 00584203208

REF 5842-32-08 **REF**

PHYSICA ZUK Unicompartmental Knee System
All-Poly Tibial Component
#1 h 8 RIGHT Medial/LEFT Lateral

2024-11-30 **LOT** 201900001 **STERILE R** 1900002



+M323005842032081/55243351900001K19R

Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

For LimaCorporate Use Only 01



(20)01(240)96171501080(10)1900001(91)1900002(92)6984



(01)08033390126588(17)241130

 Single use
  Do not use if package is damaged
  Keep away from sunlight
  temperature limit
 Use By
  **REF** Reference Number
  Keep dry
 **LOT** Lot Number
  Caution
  **STERILE EO** Sterilized using Ethylene Oxide
 IFU Consult instructions for use
  Manufacturer
  **STERILE R** Sterilized using irradiation

REF 5842-32-08 All-Poly Tibial Component
PHYSICA ZUK UHMWPE+AISI 316L
Unicompartmental Knee System

#1 h 8

RIGHT Medial/LEFT Lateral
Femoral Size A,B,C,D,E,F,G



EN All-Poly Tibial Comp DE All-Poly Tibiakomp.
 IT Comp. Tib. All-Poly FR Comp.Tibial All-Poly
 ES Comp,Tibial All-Poly PT Comp.Tibial All-Poly

Cemented Use Only

LOT 201900001 **qty:** 1 **STERILE R** Limacorporate S.p.A.

EDI	SIZE	REF
00584203210	#1 H. 10mm	5842-32-10
00584203212	#1 H. 12mm	5842-32-12
00584203214	#1 H. 14mm	5842-32-14
00584204208	#2 H. 8mm	5842-42-08
00584204210	#2 H. 10mm	5842-42-10
00584204212	#2 H. 12mm	5842-42-12
00584204214	#2 H. 14mm	5842-42-14
00584205208	#3 H. 8mm	5842-52-08
00584205210	#3 H. 10mm	5842-52-10
00584205212	#3 H. 12mm	5842-52-12
00584205214	#3 H. 14mm	5842-52-14
00584206208	#4 H. 8mm	5842-62-08
00584206210	#4 H. 10mm	5842-62-10
00584206212	#4 H. 12mm	5842-62-12
00584206214	#4 H. 14mm	5842-62-14
00584207208	#5 H. 8mm	5842-72-08
00584207210	#5 H. 10mm	5842-72-10
00584207212	#5 H. 12mm	5842-72-12
00584207214	#5 H. 14mm	5842-72-14
00584208208	#6 H. 8mm	5842-82-08
00584208210	#6 H. 10mm	5842-82-10
00584208212	#6 H. 12mm	5842-82-12
00584208214	#6 H. 14mm	5842-82-14

SECTION 2 – EXAMPLE OF LABELS FOR FEMORAL COMPONENTS (CoCrMo + PMMA)

EDI 00584201101
REF 5842-11-01
 PHYSICA ZUK Unicompartmental Knee System
 Femoral Component - High Flex
 LEFT Medial/RIGHT Lateral
 #A
 2024-11-30 LOT 201900001 STERILE R 1900002



Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

For LimaCorporate Use Only 01



(20)01(240)96110095A00(10)1900001(91)1900002(92)6543
 (01)08033390126557(17)241130

Single use Do not use if package is damaged Keep away from sunlight temperature mt
 Use By REF Reference Number Keep dry
 LOT Lot Number IFU Consult instructions for use Caution Manufacturer STERILE R Sterilized using Ethylene Oxide Sterilized using irradiation

REF 5842-11-01 Femoral Component - High Flex
 PHYSICA ZUK CoCrMo + PMMA
 Unicompartmental Knee System
 #A
 LEFT Medial/RIGHT Lateral



EN Femoral Component DE Femurkomponente
 Cemented Use Only IT Componente Femorale FR Composant Fémoral
 ES Componente Femoral PT Componente Femoral
 LOT 201900001
 2024-11-30 qty: 1 STERILE R Limacorporate S.p.A.

EDI	SIZE	REF
00584201201	#B	5842-12-01
00584201301	#C	5842-13-01
00584201401	#D	5842-14-01
00584201501	#E	5842-15-01
00584201601	#F	5842-16-01
00584201701	#G	5842-17-01

EDI

EDI 00584201102

REF 5842-11-02 **REF**

PHYSICA ZUK Unicompartamental Knee System
Femoral Component - High Flex
#A RIGHT Medial/LEFT Lateral

2024-11-30 **LOT** 201900001 **STERILE R** 1900002



*M323005842011021900001900001K191

Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy



For LimaCorporate Use Only 01



(20)01(240)96110091A00(10)1900001(91)1900002(92)6542



(01)08033390115865(17)241130

Single use	Do not use if package is damaged	Keep away from sunlight	Temperature mit
Use By	Reference Number	Keep dry	
Lot Number	Caution	Sterilized using Ethylene Oxide	
IFU Consult instructions for use	Manufacturer	Sterilized using irradiation	

REF 5842-11-02 Femoral Component - High Flex
PHYSICA ZUK CoCrMo + PMMA
Unicompartamental Knee System

#A

RIGHT Medial/LEFT Lateral

Femoral Component	Femurkomponente
Componente Femorale	Composant Fémoral
Componente Femoral	Componente Femoral

Cemented Use Only

LOT 201900001 **STERILE R** **qty:** 1 **Limacorporate S.p.A.**

EDI	SIZE	REF
00584201202	#B	5842-12-02
00584201302	#C	5842-13-02
00584201402	#D	5842-14-02
00584201502	#E	5842-15-02
00584201602	#F	5842-16-02
00584201702	#G	5842-17-02

EDI 00584200101

REF 5842-01-01

PHYSICA ZUK Unicompartmental Knee System
Tibial Component
#1 LEFT Medial/RIGHT Lateral
2024-11-30 **LOT** 201900001 **STERILE R** 1900002



+M323005842001011535243351900001K19G

Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

For LimaCorporate Use Only 01



(20)01(240)96120151010(10)1900001(91)1900002(92)6320



(01)08033390126540(17)241130

Single use Do not use if package is damaged Keep away from sunlight temperature mit
Use By REF Reference Number Keep dry
LOT Lot Number Caution **STERILE EO** Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer **STERILE R** Sterilized using irradiation

REF 5842-01-01 Tibial Component
PHYSICA ZUK Ti6Al4V + PMMA
Unicompartmental Knee System **#1**
LEFT Medial/RIGHT Lateral



EN Tibial Component **DE** Tibia-Komponente
IT Componente Tibiale **FR** Composant Tibial
ES Componente Tibial **PT** Componente Tibiall

Cemented Use Only **LOT** 201900001
2024-11-30 **qty:** 1 **STERILE R** Limacorporate S.p.A.

EDI	SIZE	REF
00584200201	#2	5842-02-01
00584200301	#3	5842-03-01
00584200401	#4	5842-04-01
00584200501	#5	5842-05-01
00584200601	#6	5842-06-01

EDI

SIZE

EDI **00584200102**
REF **5842-01-02** REF
PHYSICA ZUK Unicompartmental Knee System
Tibial Component
#1 RIGHT Medial/LEFT Lateral
2024-11-30 LOT 201900001 STERILE R 1900002

Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy 

EDI	SIZE	REF
00584200202	#2	5842-02-02
00584200302	#3	5842-03-02
00584200402	#4	5842-04-02
00584200502	#5	5842-05-02
00584200602	#6	5842-06-02

For LimaCorporate Use Only 01

(20)01(240)96120151020(10)1900001(91)1900002(92)6325

(01)08033390126564(17)241130
Single use Do not use if package is damaged Keep away from sunlight temperature mit
Use By REF Reference Number Keep dry
LOT Lot Number Caution STERILE EO Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer STERILE R Sterilized using irradiation

REF **5842-01-02** Tibial Component
PHYSICA ZUK
Unicompartmental Knee System
#1
RIGHT Medial/LEFT Lateral

EN Tibial Component DE Tibia-Komponente
Cemented Use Only IT Componente Tibiale FR Composant Tibial
LOT 201900001 ES Componente Tibial PT Componente Tibiall
2024-11-30 qty: 1 STERILE R Limacorporate S.p.A.

SECTION 4 – EXAMPLE OF LABELS FOR ARTICULAR SURFACE (UHMWPE)

EDI 00584202108  (01)08033390126571
(17)241130
(10)1900001
(91)1900002

REF 5842-21-08

PHYSICA ZUK Unicompartmental Knee System
Articular Surface
#1 h 8 Femoral Size A,B,C,D,E,F,G
2024-11-30 **LOT** 201900001 **STERILE** R 1900002



Limacorporate S.p.A. Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy **CE**
0123

For LimaCorporate Use Only 01



(20)01(240)96130501080(10)1900001(91)1900002(92)6603



(01)08033390126571(17)241130

Single use Do not use if package is damaged Keep away from sunlight temperature mit
Use By REF Reference Number Keep dry
LOT Lot Number Caution **STERILE** EO Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer **STERILE** R Sterilized using irradiation

REF 5842-21-08 Articular Surface UHMWPE
PHYSICA ZUK
Unicompartmental Knee System

#1 h 8
Femoral Size A,B,C,D,E,F,G



EN Articular Surface **DE** Gleitfläche
IT Superficie Articolare **FR** Surface Articulaire
ES Superficie Articular **PT** Surface Articulaire

LOT 201900001
2024-11-30 **qty:** 1 **STERILE** R Limacorporate S.p.A.

EDI	SIZE	REF
00584202109	#1 H. 9mm	5842-21-09
00584202110	#1 H. 10mm	5842-21-10
00584202111	#1 H. 11mm	5842-21-11
00584202112	#1 H. 12mm	5842-21-12
00584202114	#1 H. 14mm	5842-21-14
00584202208	#2 H. 8mm	5842-22-08
00584202209	#2 H. 9mm	5842-22-09
00584202210	#2 H. 10mm	5842-22-10
00584202211	#2 H. 11mm	5842-22-11
00584202212	#2 H. 12mm	5842-22-12
00584202214	#2 H. 14mm	5842-22-14
00584202308	#3 H. 8mm	5842-23-08
00584202309	#3 H. 9mm	5842-23-09
00584202310	#3 H. 10mm	5842-23-10
00584202311	#3 H. 11mm	5842-23-11
00584202312	#3 H. 12mm	5842-23-12
00584202314	#3 H. 14mm	5842-23-14
00584202408	#4 H. 8mm	5842-24-08
00584202409	#4 H. 9mm	5842-24-09
00584202410	#4 H. 10mm	5842-24-10
00584202411	#4 H. 11mm	5842-24-11
00584202412	#4 H. 12mm	5842-24-12
00584202414	#4 H. 14mm	5842-24-14
00584202508	#5 H. 8mm	5842-25-08
00584202509	#5 H. 9mm	5842-25-09
00584202510	#5 H. 10mm	5842-25-10
00584202511	#5 H. 11mm	5842-25-11
00584202512	#5 H. 12mm	5842-25-12
00584202514	#5 H. 14mm	5842-25-14
00584202608	#6 H. 8mm	5842-26-08
00584202609	#6 H. 9mm	5842-26-09
00584202610	#6 H. 10mm	5842-26-10
00584202611	#6 H. 11mm	5842-26-11
00584202612	#6 H. 12mm	5842-26-12
00584202614	#6 H. 14mm	5842-26-14