

Labels of components of Physica system

Label of the first code of each kind of device for Physica knee system components is following reported. For each label, all the information which changes in the other labels of the same device is showed.

The labels are reported as follows:

Section 1	Example of labels for KR femoral components
Section 2	Example of labels for CR femoral components
Section 3	Example of labels for PS femoral components
Section 4	Example of labels for pegs for PS femoral components
Section 5	Example of labels for tibial plates
Section 6	Example of labels for KR liners
Section 7	Example of labels for CR liners
Section 8	Example of labels for PS liners
Section 9	Example of labels for tibial stems
Section 10	Example of labels for patellar prosthesis

STERILIZATION – DATA ON THE LABELS

The information about the sterilization method for the device is reported on the labels. The external cardboard box (first label in each of the following pages) indicates only the possible sterilization methods:



- **STERILE R:** sterilization by **irradiation** (Beta rays or Gamma rays according to the applicable sterilization method).




- **STERILE EO:** sterilization by **ethylene oxide**.

Always refer to device labels for information about the sterilization method.

SECTION 1 – EXAMPLE OF LABELS FOR KR FEMORAL COMPONENTS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)



Code
REF **6511.09.110**

Size
#1

R/L
R

Sterilization
STERILE EO

Cemented Use Only
2022-07-31

PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro
Right Femoral Comp. EN
Comp. Femorale Destro IT
Comp. Femoral Derecho ES
Rechte Femurkomponente DE
Comp. Femorale Droite FR
Comp. Femoral Direito PT
Limacorporate S.p.a.

REF **6511.09.110**
PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro
2022-07-31 LOT 2017 AS098 STERILE EO 1700214



(20100)(240)96511091100(10)17AS098(91)1700214(52)1435
(01)08033390083089(17)220731

Single use, Do not use if package is damaged, Keep away from sunlight, Temperature limit, Use By, Reference Number, Keep dry, Lot Number, Caution, Sterilized using Ethylene Oxide, IFU Consult instructions for use, Manufacturer, Sterilized using irradiation.

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF **6511.09.110**
PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro
2022-07-31 LOT 2017 AS098 STERILE EO 1700214



(20100)(240)96511091100(10)17AS098(91)1700214(52)1435
(01)08033390083089(17)220731

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



WARNING

COUPLE ONLY WITH LINERS:

#1-2-3 KR RIGHT

6531.50.1XX
6531.50.2XX
6531.50.3XX

Couplings

XX:ALL THICKNESSES

Code	Size	R / L	Coupling		Applicable Sterilization method
6511.09.110	#1	R (Right)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays
6511.09.120	#2	R (Right)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6511.09.130	#3	R (Right)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6511.09.140	#4	R (Right)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6511.09.150	#5	R (Right)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6511.09.160	#6	R (Right)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6511.09.170	#7	R (Right)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6511.09.180	#8	R (Right)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6511.09.190	#9	R (Right)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6511.09.1A0	#10	R (Right)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays
6511.09.510	#1	L (Left)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays

Code	Size	R / L	Coupling		Applicable Sterilization method
6511.09.520	#2	L (Left)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6511.09.530	#3	L (Left)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6511.09.540	#4	L (Left)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6511.09.550	#5	L (Left)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6511.09.560	#6	L (Left)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6511.09.570	#7	L (Left)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6511.09.580	#8	L (Left)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6511.09.590	#9	L (Left)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6511.09.5A0	#10	L (Left)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays

Example of labels for the first code

REF 6511.09.110
PHYSICA KR
CoCrMo



#1
R

Cemented Use Only
2022-07-31 qty: 1 **STERILE EO** **Limacorporate S.p.a.**

PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro
Right Femoral Comp. **EN**
Comp. Femorale Destro **IT**
Comp. Femoral Derecho **ES**
Rechte Femurkomponente **DE**
Comp. Femorale Droite **FR**
Comp. Femoral Direito **PT**

REF 6511.09.110
PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro

2022-07-31 **LOT** 2017 AS098 **STERILE EO** 1700214



(20)00(240)96511091100(10)17AS098(91)1700214(92)1435
(01)08033390083089(17)220731

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution **STERILE EO** Sterilized using Ethylene Oxide

IFU Consult instructions for use Manufacturer **STERILE R** Sterilized using irradiation

**Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy**

CE 0123

REF 6511.09.110
PHYSICA KR - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro

2022-07-31 **LOT** 2017 AS098 **STERILE EO** 1700214



(20)00(240)96511091100(10)17AS098(91)1700214(92)1435
(01)08033390083089(17)220731

**Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy**

CE 0123



WARNING

COUPLE ONLY WITH LINERS:

#1-2-3 KR RIGHT


6531.50.1XX
6531.50.2XX
6531.50.3XX

XX:ALL THICKNESSES

SECTION 2 – EXAMPLE OF LABELS FOR CR FEMORAL COMPONENTS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)



PHYSICA CR - CoCrMo
Femoral Component / Comp. Femorale
#1 CR Cemented - Right / Destro
Right Femoral Comp.
Comp. Femorale Destro
Comp. Femoral Derecho
Rechte Femurkomponente
Comp. Femorale Droite
Comp. Femoral Direito
Limacorporate S.p.a.

Cemented Use Only
2020-10-31 qty: 1 STERILE EO
Sterilization

Single use
Do not use if package is damaged
Keep away from sunlight
Temperature limit
Use By
Reference Number
Keep dry
Lot Number
Caution
Sterilized using Ethylene Oxide
IFU Consult instructions for use
Manufacturer
Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy
CE 0123

REF 6513.09.110
PHYSICA CR - CoCrMo
Femoral Component / Comp. Femorale
#1 CR Cemented - Right / Destro
2020-10-31 LOT 2015 AS0EC STERILE EO 1500322
00

REF 6513.09.110
PHYSICA CR - CoCrMo
Femoral Component / Comp. Femorale
#1 CR Cemented - Right / Destro
2020-10-31 LOT 2015 AS0EC STERILE EO 1500322
00

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy
CE 0123


WARNING

COUPLE ONLY WITH LINERS:

#1-2-3 CR

6533.50.1XX
6533.50.2XX
6533.50.3XX
Couplings

XX: ALL THICKNESSES

Code	Size	R / L	Coupling		Applicable Sterilization method
6513.09.110	#1	R (Right)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays
6513.09.120	#2	R (Right)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6513.09.130	#3	R (Right)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6513.09.140	#4	R (Right)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6513.09.150	#5	R (Right)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6513.09.160	#6	R (Right)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6513.09.170	#7	R (Right)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6513.09.180	#8	R (Right)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6513.09.190	#9	R (Right)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6513.09.1A0	#10	R (Right)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays
6513.09.510	#1	L (Left)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays

Code	Size	R / L	Coupling		Applicable Sterilization method
6513.09.520	#2	L (Left)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6513.09.530	#3	L (Left)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6513.09.540	#4	L (Left)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6513.09.550	#5	L (Left)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6513.09.560	#6	L (Left)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6513.09.570	#7	L (Left)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6513.09.580	#8	L (Left)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6513.09.590	#9	L (Left)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6513.09.5A0	#10	L (Left)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays

Example of labels for the first code
REF 6513.09.110

 PHYSICA CR
 CoCrMo

#1
R

Cemented Use Only

2020-10-31 qty: 1

STERILE EO

 PHYSICA CR - CoCrMo
 Femoral Component / Comp. Femorale
 #1 CR Cemented - Right / Destro
Right Femoral Comp. EN
Comp. Femorale Destro IT
Comp. Femoral Derecho ES
Rechte Femurkomponente DE
Comp. Femorale Droite FR
Comp. Femoral Direito PT
Limacorporate S.p.a.
REF 6513.09.110

PHYSICA CR - CoCrMo

Femoral Component / Comp. Femorale

#1 CR Cemented - Right / Destro

2020-10-31 LOT 2015 AS0EC STERILE EO 1500322



(01)08033390087407(17)201031


 Limacorporate S.p.A. - Via Nazionale, 52
 33038 Villanova di San Daniele - Udine - Italy

CE 0123


REF 6513.09.110

PHYSICA CR - CoCrMo

Femoral Component / Comp. Femorale

#1 CR Cemented - Right / Destro

2020-10-31 LOT 2015 AS0EC STERILE EO 1500322



(01)08033390087407(17)201031

 Limacorporate S.p.A. - Via Nazionale, 52
 33038 Villanova di San Daniele - Udine - Italy 0123

WARNING
COUPLE ONLY WITH LINERS:

#1-2-3 CR

6533.50.1XX

6533.50.2XX



6533.50.3XX

XX: ALL THICKNESSES

SECTION 3 – EXAMPLE OF LABELS FOR PS FEMORAL COMPONENTS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)



Code	Size	R / L	Coupling		Applicable Sterilization method
6515.09.110	#1	R (Right)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays
6515.09.120	#2	R (Right)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6515.09.130	#3	R (Right)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6515.09.140	#4	R (Right)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6515.09.150	#5	R (Right)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6515.09.160	#6	R (Right)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6515.09.170	#7	R (Right)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6515.09.180	#8	R (Right)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6515.09.190	#9	R (Right)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6515.09.1A0	#10	R (Right)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays
6515.09.510	#1	L (Left)	#1-2-3	1XX,2XX,3XX	EtO, Gamma rays

Code	Size	R / L	Coupling		Applicable Sterilization method
6515.09.520	#2	L (Left)	#1-2-3-4	1XX,2XX,3XX,4XX	EtO, Gamma rays
6515.09.530	#3	L (Left)	#1-2-3-4-5	1XX,2XX,3XX,4XX,5XX	EtO, Gamma rays
6515.09.540	#4	L (Left)	#2-3-4-5-6	2XX,3XX,4XX,5XX,6XX	EtO, Gamma rays
6515.09.550	#5	L (Left)	#3-4-5-6-7	3XX,4XX,5XX,6XX,7XX	EtO, Gamma rays
6515.09.560	#6	L (Left)	#4-5-6-7-8	4XX,5XX,6XX,7XX,8XX	EtO, Gamma rays
6515.09.570	#7	L (Left)	#5-6-7-8-9	5XX,6XX,7XX,8XX,9XX	EtO, Gamma rays
6515.09.580	#8	L (Left)	#6-7-8-9-A	6XX,7XX,8XX,9XX,AXX	EtO, Gamma rays
6515.09.590	#9	L (Left)	#7-8-9-A	7XX,8XX,9XX,AXX	EtO, Gamma rays
6515.09.5A0	#10	L (Left)	#8-9-A	8XX,9XX,AXX	EtO, Gamma rays

Example of labels for the first code

REF 6515.09.110
PHYSICA PS - CoCrMo
CoCrMo



#1
R

Cemented Use Only
2022-03-31 qty: 1 STERILE EO Limacorporate S.p.a.

PHYSICA PS - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro
Right Femoral Comp. **EN**
Comp. Femorale Destro **IT**
Comp. Femoral Derecho **ES**
Rechte Femurkomponente **DE**
Comp. Femorale Droite **FR**
Comp. Femoral Direito **PT**

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution Sterilized using Ethylene Oxide

IFU Consult instructions for use Manufacturer Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF 6515.09.110
PHYSICA PS - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro

2022-03-31 LOT 2017 AS03R STERILE EO 1700081



(20)00(240)96515091100(10)17AS03R(91)1700081(92)1437
(01)08033390095495(17)220331

REF 6515.09.110
PHYSICA PS - CoCrMo
Femoral Component / Comp. Femorale
#1Cemented - Right / Destro

2022-03-31 LOT 2017 AS03R STERILE EO 1700081



(20)00(240)96515091100(10)17AS03R(91)1700081(92)1437
(01)08033390095495(17)220331

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123


WARNING

COUPLE ONLY 361
WITH LINERS:

#1-2-3 PS

6535.50.1XX

6535.50.2XX

6535.50.3XX

XX:ALL THICKNESSES

SECTION 4 – EXAMPLE OF LABELS FOR PEGS FOR PS FEMORAL COMPONENTS

REF 6515.09.900

PHYSICA PS
CoCrMo



2 Pieces in 1 Pack
2021-11-30 **qty:** 1 **STERILE R**

PHYSICA PS - CoCrMo
PS Femoral Peg
Perno Femorale PS
PS Femoral Peg EN
Perno Femorale PS IT
Tetón Del Femoral PS ES
PS Femuzapfen DE
Plot Fémoral PS FR
Peg Femorais PS PT

Limacorporate S.p.a.

REF 6515.09.900

PHYSICA PS - CoCrMo
PS Femoral Peg
Perno Femorale PS

2021-11-30 **LOT** 2016 15127 **STERILE R** 1600296



(20)00(240)96515099000(10)1615127(91)1600296(92)1438



(01)08033390100021(17)211130



 Single use

 Do not use if package is damaged

 Keep away from sunlight

 Temperature limit

 Use By

 **REF** Reference Number

 Keep dry

 **LOT** Lot Number

 **CAUTION** Caution

 **STERILE EO** Sterilized using Ethylene Oxide

 IFU Consult instructions for use

 **Manufacturer**

 **STERILE R** Sterilized using irradiation


Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy


0123



REF 6515.09.900

PHYSICA PS - CoCrMo
PS Femoral Peg
Perno Femorale PS

2021-11-30 **LOT** 2016 15127 **STERILE R** 1600296



(20)00(240)96515099000(10)1615127(91)1600296(92)1438



(01)08033390100021(17)211130




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy


0123


No warnings

Code	Applicable Sterilization method
6515.09.900	Beta rays, Gamma rays

SECTION 5 – EXAMPLE OF LABELS FOR TIBIAL PLATES

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)



No warnings

Code	Size	Applicable Sterilization method
6522.15.010	#1	Beta rays
6522.15.020	#2	Beta rays
6522.15.030	#3	Beta rays
6522.15.040	#4	Beta rays
6522.15.050	#5	Beta rays
6522.15.060	#6	Beta rays
6522.15.070	#7	Beta rays
6522.15.080	#8	Beta rays
6522.15.090	#9	Beta rays
6522.15.0A0	#10	Beta rays

Example of labels for the first code

REF 6522.15.010

PHYSICA
Ti6Al4V + UHMWPE



#1

Cemented Use Only

2022-08-31 qty: 1 STERILE R Limacorporate S.p.a.

PHYSICA Knee - Ti6Al4V + UHMWPE
Fixed Tibial Plate/Piatto Tib.Fisso
#1 Cemented / Cementato
Fixed Tibial Plate **EN**
Piatto Tibiale Fisso **IT**
Plato Tibial Fijo **ES**
Fix Tibiplateau **DE**
Plateau Tibial Fixe **FR**
Prato Tibial Fijo **PT**

REF 6522.15.010

PHYSICA Knee - Ti6Al4V + UHMWPE
Fixed Tibial Plate/Piatto Tib.Fisso

#1 Cemented / Cementato

2022-08-31 LOT 2016 07574 STERILE R 1700256



(20)00(240)96522150100(10)1607574(91)1700256(92)1439

(01)08033390088350(17)220831



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 6522.15.010

PHYSICA Knee - Ti6Al4V + UHMWPE
Fixed Tibial Plate/Piatto Tib.Fisso

#1 Cemented / Cementato

2022-08-31 LOT 2016 07574 STERILE R 1700256



(20)00(240)96522150100(10)1607574(91)1700256(92)1439

(01)08033390088350(17)220831



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy


CE
0123

No warnings

SECTION 6 – EXAMPLE OF LABELS FOR KR LINERS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)



Code 6531.50.110
PHYSICA KR - UHMWPE

Size #1
R/L R
H 10 H

2022-06-30 qty: 1 **STERILE** EO
Sterilization

Limacorporate S.p.a.

REF 6531.50.110
PHYSICA KR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm - Right/Destro
Right Tibial Liner **EN**
Inserto Tibiale Destro **IT**
Inserto Tibial Derecho **ES**
Rechtes Tibial Inlay **DE**
Insert Tibial Droite **FR**
Polietileno Tibia Direito **PT**

2022-06-30 **LOT** 2016 AT1KH **STERILE** EO 1700207

(20)00(240)96531501100(10)16AT1KH(91)1700207(92)1440
(01)08033390087629(17)220630

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF 6531.50.110
PHYSICA KR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm - Right/Destro
2022-06-30 **LOT** 2016 AT1KH **STERILE** EO 1700207

(20)00(240)96531501100(10)16AT1KH(91)1700207(92)1440
(01)08033390087629(17)220630

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



WARNING

COUPLE ONLY ³¹¹
WITH FEMORAL COMPONENTS:
#1-2-3 KR RIGHT

6511.09.110
6511.09.120
6511.09.130

Couplings

Code	Size	R / L	H	Coupling		Applicable Sterilization method
6531.50.110	#1	R (Right)	10	#1-2-3	110,120,130	EtO
6531.50.111	#1	R (Right)	11	#1-2-3	110,120,130	EtO
6531.50.112	#1	R (Right)	12	#1-2-3	110,120,130	EtO
6531.50.114	#1	R (Right)	14	#1-2-3	110,120,130	EtO
6531.50.116	#1	R (Right)	16	#1-2-3	110,120,130	EtO
6531.50.120	#1	R (Right)	20	#1-2-3	110,120,130	EtO
6531.50.210	#2	R (Right)	10	#1-2-3-4	110,120,130,140	EtO
6531.50.211	#2	R (Right)	11	#1-2-3-4	110,120,130,140	EtO
6531.50.212	#2	R (Right)	12	#1-2-3-4	110,120,130,140	EtO
6531.50.214	#2	R (Right)	14	#1-2-3-4	110,120,130,140	EtO

Code	Size	R / L	H	Coupling		Applicable Sterilization method
6531.50.216	#2	R (Right)	16	#1-2-3-4	110,120,130,140	EtO
6531.50.220	#2	R (Right)	20	#1-2-3-4	110,120,130,140	EtO
6531.50.310	#3	R (Right)	10	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.311	#3	R (Right)	11	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.312	#3	R (Right)	12	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.314	#3	R (Right)	14	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.316	#3	R (Right)	16	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.320	#3	R (Right)	20	#1-2-3-4-5	110,120,130,140,150	EtO
6531.50.410	#4	R (Right)	10	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.411	#4	R (Right)	11	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.412	#4	R (Right)	12	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.414	#4	R (Right)	14	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.416	#4	R (Right)	16	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.420	#4	R (Right)	20	#2-3-4-5-6	120,130,140,150,160	EtO
6531.50.510	#5	R (Right)	10	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.511	#5	R (Right)	11	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.512	#5	R (Right)	12	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.514	#5	R (Right)	14	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.516	#5	R (Right)	16	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.520	#5	R (Right)	20	#3-4-5-6-7	130,140,150,160,170	EtO
6531.50.610	#6	R (Right)	10	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.611	#6	R (Right)	11	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.612	#6	R (Right)	12	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.614	#6	R (Right)	14	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.616	#6	R (Right)	16	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.620	#6	R (Right)	20	#4-5-6-7-8	140,150,160,170,180	EtO
6531.50.710	#7	R (Right)	10	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.711	#7	R (Right)	11	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.712	#7	R (Right)	12	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.714	#7	R (Right)	14	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.716	#7	R (Right)	16	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.720	#7	R (Right)	20	#5-6-7-8-9	150,160,170,180,190	EtO
6531.50.810	#8	R (Right)	10	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.811	#8	R (Right)	11	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.812	#8	R (Right)	12	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.814	#8	R (Right)	14	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.816	#8	R (Right)	16	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.820	#8	R (Right)	20	#6-7-8-9-10	160,170,180,190,1A0	EtO
6531.50.910	#9	R (Right)	10	#7-8-9-10	170,180,190,1A0	EtO
6531.50.911	#9	R (Right)	11	#7-8-9-10	170,180,190,1A0	EtO
6531.50.912	#9	R (Right)	12	#7-8-9-10	170,180,190,1A0	EtO

Code	Size	R / L	H	Coupling		Applicable Sterilization method
6531.50.914	#9	R (Right)	14	#7-8-9-10	170,180,190,1A0	EtO
6531.50.916	#9	R (Right)	16	#7-8-9-10	170,180,190,1A0	EtO
6531.50.920	#9	R (Right)	20	#7-8-9-10	170,180,190,1A0	EtO
6531.50.A10	#10	R (Right)	10	#8-9-10	180,190,1A0	EtO
6531.50.A11	#10	R (Right)	11	#8-9-10	180,190,1A0	EtO
6531.50.A12	#10	R (Right)	12	#8-9-10	180,190,1A0	EtO
6531.50.A14	#10	R (Right)	14	#8-9-10	180,190,1A0	EtO
6531.50.A16	#10	R (Right)	16	#8-9-10	180,190,1A0	EtO
6531.50.A20	#10	R (Right)	20	#8-9-10	180,190,1A0	EtO
6532.50.110	#1	L(Left)	10	#1-2-3	510,520,530	EtO
6532.50.111	#1	L(Left)	11	#1-2-3	510,520,530	EtO
6532.50.112	#1	L(Left)	12	#1-2-3	510,520,530	EtO
6532.50.114	#1	L(Left)	14	#1-2-3	510,520,530	EtO
6532.50.116	#1	L(Left)	16	#1-2-3	510,520,530	EtO
6532.50.120	#1	L(Left)	20	#1-2-3	510,520,530	EtO
6532.50.210	#2	L(Left)	10	#1-2-3-4	510,520,530,540	EtO
6532.50.211	#2	L(Left)	11	#1-2-3-4	510,520,530,540	EtO
6532.50.212	#2	L(Left)	12	#1-2-3-4	510,520,530,540	EtO
6532.50.214	#2	L(Left)	14	#1-2-3-4	510,520,530,540	EtO
6532.50.216	#2	L(Left)	16	#1-2-3-4	510,520,530,540	EtO
6532.50.220	#2	L(Left)	20	#1-2-3-4	510,520,530,540	EtO
6532.50.310	#3	L(Left)	10	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.311	#3	L(Left)	11	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.312	#3	L(Left)	12	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.314	#3	L(Left)	14	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.316	#3	L(Left)	16	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.320	#3	L(Left)	20	#1-2-3-4-5	510,520,530,540,550	EtO
6532.50.410	#4	L(Left)	10	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.411	#4	L(Left)	11	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.412	#4	L(Left)	12	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.414	#4	L(Left)	14	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.416	#4	L(Left)	16	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.420	#4	L(Left)	20	#2-3-4-5-6	520,530,540,550,560	EtO
6532.50.510	#5	L(Left)	10	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.511	#5	L(Left)	11	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.512	#5	L(Left)	12	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.514	#5	L(Left)	14	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.516	#5	L(Left)	16	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.520	#5	L(Left)	20	#3-4-5-6-7	530,540,550,560,570	EtO
6532.50.610	#6	L(Left)	10	#4-5-6-7-8	540,550,560,570,580	EtO
6532.50.611	#6	L(Left)	11	#4-5-6-7-8	540,550,560,570,580	EtO

Code	Size	R / L	H	Coupling		Applicable Sterilization method
6532.50.612	#6	L(Left)	12	#4-5-6-7-8	540,550,560,570,580	EtO
6532.50.614	#6	L(Left)	14	#4-5-6-7-8	540,550,560,570,580	EtO
6532.50.616	#6	L(Left)	16	#4-5-6-7-8	540,550,560,570,580	EtO
6532.50.620	#6	L(Left)	20	#4-5-6-7-8	540,550,560,570,580	EtO
6532.50.710	#7	L(Left)	10	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.711	#7	L(Left)	11	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.712	#7	L(Left)	12	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.714	#7	L(Left)	14	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.716	#7	L(Left)	16	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.720	#7	L(Left)	20	#5-6-7-8-9	550,560,570,580,590	EtO
6532.50.810	#8	L(Left)	10	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.811	#8	L(Left)	11	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.812	#8	L(Left)	12	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.814	#8	L(Left)	14	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.816	#8	L(Left)	16	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.820	#8	L(Left)	20	#6-7-8-9-10	560,570,580,590,5A0	EtO
6532.50.910	#9	L(Left)	10	#7-8-9-10	570,580,590,5A0	EtO
6532.50.911	#9	L(Left)	11	#7-8-9-10	570,580,590,5A0	EtO
6532.50.912	#9	L(Left)	12	#7-8-9-10	570,580,590,5A0	EtO
6532.50.914	#9	L(Left)	14	#7-8-9-10	570,580,590,5A0	EtO
6532.50.916	#9	L(Left)	16	#7-8-9-10	570,580,590,5A0	EtO
6532.50.920	#9	L(Left)	20	#7-8-9-10	570,580,590,5A0	EtO
6532.50.A10	#10	L(Left)	10	#8-9-10	580,590,5A0	EtO
6532.50.A11	#10	L(Left)	11	#8-9-10	580,590,5A0	EtO
6532.50.A12	#10	L(Left)	12	#8-9-10	580,590,5A0	EtO
6532.50.A14	#10	L(Left)	14	#8-9-10	580,590,5A0	EtO
6532.50.A16	#10	L(Left)	16	#8-9-10	580,590,5A0	EtO
6532.50.A20	#10	L(Left)	20	#8-9-10	580,590,5A0	EtO

Example of labels for the first code

REF 6531.50.110
PHYSICA KR
UHMWPE



#1
R
H 10

PHYSICA KR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm - Right/Destro
Right Tibial Liner **EN**
Inserto Tibiale Destro **IT**
Inserto Tibial Derecho **ES**
Rechtes Tibial Inlay **DE**
Insert Tibial Droite **FR**
Polietileno Tibia Direito **PT**
Limacorporate S.p.a.

2022-06-30 **qty:** 1 **STERILE** **EO**

REF 6531.50.110
PHYSICA KR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm - Right/Destro

2022-06-30 **LOT** 2016 AT1KH **STERILE** **EO** 1700207



(20)00(240)96531501100(10)16AT1KH(91)1700207(92)1440
(01)08033390087629(17)220630

 Single use	 Do not use if package is damaged	 Keep away from sunlight	 Temperature limit
 Use By	REF Reference Number	 Keep dry	
LOT Lot Number	 Caution	STERILE EO Sterilized using Ethylene Oxide	
 IFU Consult instructions for use	 Manufacturer	STERILE R Sterilized using irradiation	

 Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123 

REF 6531.50.110
PHYSICA KR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm - Right/Destro

2022-06-30 **LOT** 2016 AT1KH **STERILE** **EO** 1700207



(20)00(240)96531501100(10)16AT1KH(91)1700207(92)1440
(01)08033390087629(17)220630

 Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy **CE** 0123

 WARNING	COUPLE ONLY 311 WITH FEMORAL COMPONENTS: #1-2-3 KR RIGHT 6511.09.110 6511.09.120 6511.09.130
---	---

SECTION 7 – EXAMPLE OF LABELS FOR CR LINERS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)

Code
REF **6533.50.110**

PHYSICA CR - UHMWPE



Size
#1

H 10 H

2020-06-30 qty. 1 **STERILE EO** Sterilization

Limacorporate S.p.a.

PHYSICA CR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

Tibial Liner **EN**
Inserto Tibiale **IT**
Inserto Tibiale **ES**
Tibial Inlay **DE**
Insert Tibial **FR**
Polietileno Tibiale **PT**

2020-06-30 **LOT** 2015 AT02Z **STERILE EO** 1500206

(20)00(240)96533501100(10)15AT02Z(91)1500206(92)1442

(01)08033390088374(17)200630

Single use, Do not use if package is damaged, Keep away from sunlight, Temperature limit, Use By, Reference Number, Keep dry, Lot Number, Caution, Sterilized using Ethylene Oxide, IFU Consult instructions for use, Manufacturer, Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF **6533.50.110**

PHYSICA CR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2020-06-30 **LOT** 2015 AT02Z **STERILE EO** 1500206

(20)00(240)96533501100(10)15AT02Z(91)1500206(92)1442

(01)08033390088374(17)200630

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



WARNING

COUPLE ONLY ³⁵¹

WITH FEMORAL COMPONENTS:
LEFT AND RIGHT

#1-2-3 CR

6513.09.110-510
6513.09.120-520
6513.09.130-530
Couplings

Code	Size	H	Coupling		Applicable Sterilization method
6533.50.110	#1	10	#1-2-3	110,120,130,510,520,530	EtO
6533.50.111	#1	11	#1-2-3	110,120,130,510,520,530	EtO
6533.50.112	#1	12	#1-2-3	110,120,130,510,520,530	EtO
6533.50.114	#1	14	#1-2-3	110,120,130,510,520,530	EtO
6533.50.116	#1	16	#1-2-3	110,120,130,510,520,530	EtO
6533.50.120	#1	20	#1-2-3	110,120,130,510,520,530	EtO
6533.50.210	#2	10	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6533.50.211	#2	11	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6533.50.212	#2	12	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6533.50.214	#2	14	#1-2-3-4	110,120,130,140,510,520,530,540	EtO

Code	Size	H	Coupling	Applicable Sterilization method
6533.50.216	#2	16	#1-2-3-4 110,120,130,140,510,520,530,540	EtO
6533.50.220	#2	20	#1-2-3-4 110,120,130,140,510,520,530,540	EtO
6533.50.310	#3	10	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.311	#3	11	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.312	#3	12	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.314	#3	14	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.316	#3	16	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.320	#3	20	#1-2-3-4-5 110,120,130,140,150,510,520,530,540,550	EtO
6533.50.410	#4	10	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.411	#4	11	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.412	#4	12	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.414	#4	14	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.416	#4	16	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.420	#4	20	#2-3-4-5-6 120,130,140,150,160,520,530,540,550,560	EtO
6533.50.510	#5	10	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.511	#5	11	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.512	#5	12	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.514	#5	14	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.516	#5	16	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.520	#5	20	#3-4-5-6-7 130,140,150,160,170,530,540,550,560,570	EtO
6533.50.610	#6	10	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.611	#6	11	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.612	#6	12	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.614	#6	14	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.616	#6	16	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.620	#6	20	#4-5-6-7-8 140,150,160,170,180,540,550,560,570,580	EtO
6533.50.710	#7	10	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.711	#7	11	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.712	#7	12	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.714	#7	14	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.716	#7	16	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.720	#7	20	#5-6-7-8-9 150,160,170,180,190,550,560,570,580,590	EtO
6533.50.810	#8	10	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.811	#8	11	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.812	#8	12	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.814	#8	14	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.816	#8	16	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.820	#8	20	#6-7-8-9-10 160,170,180,190,1A0,560,570,580,590,5A0	EtO
6533.50.910	#9	10	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO
6533.50.911	#9	11	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO

Code	Size	H	Coupling		Applicable Sterilization method
6533.50.912	#9	12	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO
6533.50.914	#9	14	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO
6533.50.916	#9	16	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO
6533.50.920	#9	20	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO
6533.50.A10	#10	10	#8-9-10	180,190,1A0,580,590,5A0	EtO
6533.50.A11	#10	11	#8-9-10	180,190,1A0,580,590,5A0	EtO
6533.50.A12	#10	12	#8-9-10	180,190,1A0,580,590,5A0	EtO
6533.50.A14	#10	14	#8-9-10	180,190,1A0,580,590,5A0	EtO
6533.50.A16	#10	16	#8-9-10	180,190,1A0,580,590,5A0	EtO
6533.50.A20	#10	20	#8-9-10	180,190,1A0,580,590,5A0	EtO

Example of labels for the first code

REF 6533.50.110
PHYSICA CR
UHMWPE



#1
H 10

PHYSICA CR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm
Tibial Liner **EN**
Inserto Tibiale **IT**
Inserto Tibiale **ES**
Tibial Inlay **DE**
Insert Tibial **FR**
Polietileno Tibiale **PT**

2020-06-30 **qty:** 1 **STERILE EO** **Limacorporate S.p.a.**

REF 6533.50.110
PHYSICA CR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2020-06-30 **LOT** 2015 AT02Z **STERILE EO** 1500206



(20)00(240)96533501100(10)15AT02Z(91)1500206(92)1442
(01)08033390088374(17)200630

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution **STERILE EO** Sterilized using Ethylene Oxide

IFU Consult instructions for use Manufacturer **STERILE** Sterilized using irradiation

**Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy**

CE 0123

REF 6533.50.110
PHYSICA CR - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2020-06-30 **LOT** 2015 AT02Z **STERILE EO** 1500206



(20)00(240)96533501100(10)15AT02Z(91)1500206(92)1442
(01)08033390088374(17)200630

**Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy**

CE 0123



WARNING

COUPLE ONLY **351**

**WITH FEMORAL COMPONENTS:
LEFT AND RIGHT**

#1-2-3 CR

6513.09.110-510
6513.09.120-520
6513.09.130-530

SECTION 8 – EXAMPLE OF LABELS FOR PS LINERS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)

Code
REF **6535.50.110**

PHYSICA PS - UHMWPE
Tibial Liner / Inserto Tibiale

Size
#1

H 10

2022-07-31 qty: 1 **STERILE EO** Sterilization
Limacorporate S.p.a.

Single use Do not use if package is damaged Keep away from sunlight Temperature limit
Use By REF Reference Number Keep dry
LOT Lot Number Caution Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF **6535.50.110**

PHYSICA PS - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2022-07-31 LOT 2016 AT1YU **STERILE EO** 1700232

(20)00(240)96535501100(10)16AT1YU(91)1700232(92)1443
(01)08033390095693(17)220731

REF **6535.50.110**

PHYSICA PS - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2022-07-31 LOT 2016 AT1YU **STERILE EO** 1700232

(20)00(240)96535501100(10)16AT1YU(91)1700232(92)1443
(01)08033390095693(17)220731

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF **6535.50.110**

PHYSICA PS - UHMWPE
Tibial Liner / Inserto Tibiale
#1 H.10mm

2022-07-31 LOT 2016 AT1YU **STERILE EO** 1700232

(20)00(240)96535501100(10)16AT1YU(91)1700232(92)1443
(01)08033390095693(17)220731

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



WARNING

COUPLE ONLY **371**

WITH FEMORAL COMPONENTS:

LEFT AND RIGHT
#1-2-3 PS

6515.09.110-510
6515.09.120-520
~~6515.09.130-530~~

Code	Size	H	Coupling		Applicable Sterilization method
6535.50.110	#1	10	#1-2-3	110,120,130,510,520,530	EtO
6535.50.111	#1	11	#1-2-3	110,120,130,510,520,530	EtO
6535.50.112	#1	12	#1-2-3	110,120,130,510,520,530	EtO
6535.50.114	#1	14	#1-2-3	110,120,130,510,520,530	EtO
6535.50.116	#1	16	#1-2-3	110,120,130,510,520,530	EtO
6535.50.120	#1	20	#1-2-3	110,120,130,510,520,530	EtO
6535.50.210	#2	10	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6535.50.211	#2	11	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6535.50.212	#2	12	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6535.50.214	#2	14	#1-2-3-4	110,120,130,140,510,520,530,540	EtO

Code	Size	H	Coupling		Applicable Sterilization method
6535.50.216	#2	16	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6535.50.220	#2	20	#1-2-3-4	110,120,130,140,510,520,530,540	EtO
6535.50.310	#3	10	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.311	#3	11	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.312	#3	12	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.314	#3	14	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.316	#3	16	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.320	#3	20	#1-2-3-4-5	110,120,130,140,150,510,520,530,540,550	EtO
6535.50.410	#4	10	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.411	#4	11	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.412	#4	12	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.414	#4	14	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.416	#4	16	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.420	#4	20	#2-3-4-5-6	120,130,140,150,160,520,530,540,550,560	EtO
6535.50.510	#5	10	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.511	#5	11	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.512	#5	12	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.514	#5	14	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.516	#5	16	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.520	#5	20	#3-4-5-6-7	130,140,150,160,170,530,540,550,560,570	EtO
6535.50.610	#6	10	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.611	#6	11	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.612	#6	12	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.614	#6	14	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.616	#6	16	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.620	#6	20	#4-5-6-7-8	140,150,160,170,180,540,550,560,570,580	EtO
6535.50.710	#7	10	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.711	#7	11	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.712	#7	12	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.714	#7	14	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.716	#7	16	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.720	#7	20	#5-6-7-8-9	150,160,170,180,190,550,560,570,580,590	EtO
6535.50.810	#8	10	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.811	#8	11	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.812	#8	12	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.814	#8	14	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.816	#8	16	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.820	#8	20	#6-7-8-9-10	160,170,180,190,1A0,560,570,580,590,5A0	EtO
6535.50.910	#9	10	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO
6535.50.911	#9	11	#7-8-9-10	170,180,190,1A0,570,580,590,5A0	EtO

Code	Size	H	Coupling	Applicable Sterilization method
6535.50.912	#9	12	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO
6535.50.914	#9	14	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO
6535.50.916	#9	16	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO
6535.50.920	#9	20	#7-8-9-10 170,180,190,1A0,570,580,590,5A0	EtO
6535.50.A10	#10	10	#8-9-10 180,190,1A0,580,590,5A0	EtO
6535.50.A11	#10	11	#8-9-10 180,190,1A0,580,590,5A0	EtO
6535.50.A12	#10	12	#8-9-10 180,190,1A0,580,590,5A0	EtO
6535.50.A14	#10	14	#8-9-10 180,190,1A0,580,590,5A0	EtO
6535.50.A16	#10	16	#8-9-10 180,190,1A0,580,590,5A0	EtO
6535.50.A20	#10	20	#8-9-10 180,190,1A0,580,590,5A0	EtO

Example of labels for the first code

REF 6535.50.110
 PHYSICA PS - UHMWPE
 Tibial Liner / Inserto Tibiale
 #1 H.10mm
 Tibial Liner **EN**
 Inserto Tibiale **IT**
 Inserto Tibiale **ES**
 Tibial Inlay **DE**
 Insert Tibial **FR**
 Polietileno Tibial **PT**



#1
H 10

2022-07-31 qty: 1 **STERILE** **EO** Limacorporate S.p.a.

REF 6535.50.110
 PHYSICA PS - UHMWPE
 Tibial Liner / Inserto Tibiale
 #1 H.10mm
 2022-07-31 **LOT** 2016 AT1YU **STERILE** **EO** 1700232



(20)00(240)96535501100(10)16AT1YU(91)1700232(92)1443
 (01)08033390095693(17)220731

 Single use
  Do not use if package is damaged
  Keep away from sunlight
  Temperature limit

 Use By
  REF Reference Number
  Keep dry

 **LOT** Lot Number
  Caution
  Sterilized using Ethylene Oxide

 IFU Consult instructions for use
  Manufacturer
  Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
 33038 Villanova di San Daniele - Udine - Italy
  0123

REF 6535.50.110
 PHYSICA PS - UHMWPE
 Tibial Liner / Inserto Tibiale
 #1 H.10mm
 2022-07-31 **LOT** 2016 AT1YU **STERILE** **EO** 1700232



(20)00(240)96535501100(10)16AT1YU(91)1700232(92)1443
 (01)08033390095693(17)220731

Limacorporate S.p.A. - Via Nazionale, 52
 33038 Villanova di San Daniele - Udine - Italy
  0123



WARNING

COUPLE ONLY **371**

WITH FEMORAL COMPONENTS:
LEFT AND RIGHT
#1-2-3 PS

6515.09.110-510
6515.09.120-520
6515.09.130-530

SECTION 9 – EXAMPLE OF LABELS FOR TIBIAL STEMS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)





No warnings

Code	Size	Applicable Sterilization method
6590.15.020	20	Beta rays, Gamma rays
6590.15.040	40	Beta rays, Gamma rays
6590.15.060	60	Beta rays, Gamma rays

Example of labels for the first code

REF **6590.15.020**
PHYSICA
Ti6Al4V



L 20

PHYSICA Knee - Ti6Al4V
Tibial Stem / Stelo Tibiale
L.20mm

Tibial Stem **EN**
Stelo Tibiale **IT**
Vastago Tibiale **ES**
Tibiaschaft **DE**
Tige Tibial **FR**
Haste Tibial **PT**

2022-02-28 qty: 1 STERILE R Limacorporate S.p.a.

Single use	Do not use if package is damaged	Keep away from sunlight	Temperature limit
Use By	REF Reference Number	Keep dry	
LOT Lot Number	Caution	STERILE EO Sterilized using Ethylene Oxide	
If U Consult instructions for use	Manufacturer	STERILE R Sterilized using irradiation	

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



REF **6590.15.020**
PHYSICA Knee - Ti6Al4V
Tibial Stem / Stelo Tibiale
L.20mm

2022-02-28 LOT 2017 00100 STERILE R 1700039



(20)00(240)96590150200(10)1700100(51)1700039(52)1444



(01)08033390088398(17)220228

REF **6590.15.020**
PHYSICA Knee - Ti6Al4V
Tibial Stem / Stelo Tibiale
L.20mm

2022-02-28 LOT 2017 00100 STERILE R 1700039



(20)00(240)96590150200(10)1700100(51)1700039(52)1444



(01)08033390088398(17)220228

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



No warnings

SECTION 10 – EXAMPLE OF LABELS FOR PATELLAR PROSTHESIS

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below)

Code
REF 6595.50.026

PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata

Size
D 26
H 7 H

Cemented Use Only
2022-02-28 qty: 1 STERILE EO Limacorporate S.p.a.

Sterilization

Single use Do not use if package is damaged Keep away from sunlight Temperature limit
Use By REF Reference Number Keep dry
LOT Lot Number Caution STERILE EO Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer STERILE EO Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy 0123

REF 6595.50.026
PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata
2022-02-28 LOT 2016 AT1NW STERILE EO 1700036

(20)00(240)96595500260(10)16AT1NW(91)1700036(92)1445
(01)08033390088404(17)220228

REF 6595.50.026
PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata
2022-02-28 LOT 2016 AT1NW STERILE EO 1700036

(20)00(240)96595500260(10)16AT1NW(91)1700036(92)1445
(01)08033390088404(17)220228

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy 0123

No warnings

Code	Size	H	Applicable Sterilization method
6595.50.026	26	7	EtO
6595.50.029	29	7	EtO
6595.50.032	32	8.5	EtO
6595.50.035	35	8.5	EtO
6595.50.038	38	10	EtO
6595.50.041	41	10	EtO

Example of labels for the first code

REF 6595.50.026
PHYSICA
UHMWPE



D 26
H 7

Cemented Use Only
2022-02-28 **qty:** 1 **STERILE EO** **Limacorporate S.p.a.**

PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata
Patellar Prosthesis EN
Protesi Rotulea IT
Rotula Protetica ES
Patella Prothese DE
Prothese de Rotule FR
Rotula Protetica PT

REF 6595.50.026
PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata
2022-02-28 **LOT** 2016 AT1NW **STERILE EO** 1700036



(20)00(240)96595500260(10)16AT1NW(91)1700036(92)1445
(01)08033390088404(17)220228

 Single use	 Do not use if package is damaged	 Keep away from sunlight	 Temperature limit
 Use By	REF Reference Number	 Keep dry	
LOT Lot Number	 Caution	STERILE EO Sterilized using Ethylene Oxide	
 IFU Consult instructions for use	 Manufacturer	STERILE R Sterilized using irradiation	

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 6595.50.026
PHYSICA Knee - UHMWPE
Patellar Prostheses/Protesi Rotulea
D.26mm H.7mm - Cemented/Cementata
2022-02-28 **LOT** 2016 AT1NW **STERILE EO** 1700036



(20)00(240)96595500260(10)16AT1NW(91)1700036(92)1445
(01)08033390088404(17)220228

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

No warnings