

## **EXAMPLE OF LABELS OF MODULUS & REVISION SYSTEMS**

Label of the first code of each kind of device for Modulus & Revision Systems is following reported. For each label, all the information which changes in the other labels of the same device is showed.

The labels are reported as follows:

Section 1	Example of labels for Revision Modular Stems
Section 2	Example of labels for Revision Necks with Locking Screw
Section 3	Example of labels for Revision LR Stems
Section 4	Example of labels for Revision LR Modular Augments
Section 5	Example of labels for Modulus stems
Section 6	Example of labels for Modulus-R stems
Section 7	Example of labels for Modulus Necks with Locking Screw

## **STERILIZATION – DATA ON THE LABELS**

The information about the sterilization method for the device is reported on the labels. The external cardboard box (first label in each of the following pages) indicates only the possible sterilization methods:



- **STERILE R:** sterilization by **irradiation** (Beta rays or Gamma rays according to the applicable sterilization method).




- **STERILE EO:** sterilization by **ethylene oxide**.

Refer to device labels for information about the sterilization method.

For the information on applicable sterilization method for each device refer to Design Dossier.

**SECTION 1 – EXAMPLE OF LABELS FOR REVISION MODULAR STEMS**
**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).



**Code** REF 3810.15.010

D=14mm L=140mm  
 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm

**Diameter** **Length**

Uncemented Stem **EN**  
 Stelo non Cementato **IT**  
 Vastago no Cementado **ES**  
 Zementlos Schaft **DE**  
 Tige non Cimentée **FR**  
 Haste não Cimentada **PT**

**Sterilization** 2022-08-31 QTY: 1 STERILE R Limacorporate S.p.a.

REF 3810.15.010 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm  
 2022-08-31 LOT 2017 09554 STERILE R 1700266

REF 3810.15.010 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm  
 2022-08-31 LOT 2017 09554 STERILE R 1700266

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy

**WARNING**

**ENG:** (test according to ISO 7206-4)  
 Do not couple to lateralized necks

**ITA:** (test secondo ISO 7206-4)  
 Non accoppiabile a colli lateralizzanti

**DEU:** (Prüfung gemäß ISO7206-4)  
 Nicht geeignet für lateralisierte Halsteile

**FRA:** (test d'après ISO 7206-4)  
 Ne pas assembler aux cols latéralisants


**SPA:** (testado según ISO7206-4)  
 No acoplable a los cuellos lateralizados

**POR:** (testado segundo ISO7206-4)  
 No Adaptar a côlos lateralizados

Code	Length	Diameter
3810.15.010	L=140 mm	D=14 mm
3810.15.020	L=200 mm	D=14 mm
3811.15.010	L=140 mm	D=15 mm
3811.15.020	L=200 mm	D=15 mm
3812.15.010	L=140 mm	D=16 mm
3812.15.020	L=200 mm	D=16 mm
3813.15.010	L=140 mm	D=17 mm
3813.15.020	L=200 mm	D=17 mm
3814.15.010	L=140 mm	D=18 mm
3814.15.020	L=200 mm	D=18 mm
3815.15.010	L=140 mm	D=19 mm

Code	Length	Diameter
3815.15.020	L=200 mm	D=19 mm
3816.15.010	L=140 mm	D=20 mm
3816.15.020	L=200 mm	D=20 mm
3817.15.010	L=140 mm	D=21 mm
3817.15.020	L=200 mm	D=21 mm
3818.15.010	L=140 mm	D=22 mm
3818.15.020	L=200 mm	D=22 mm
3819.15.010	L=140 mm	D=23 mm
3819.15.020	L=200 mm	D=23 mm
3820.15.010	L=140 mm	D=24 mm
3820.15.020	L=200 mm	D=24 mm
3821.15.010	L=140 mm	D=25 mm
3821.15.020	L=200 mm	D=25 mm
3822.15.010	L=140 mm	D=26 mm
3822.15.020	L=200 mm	D=26 mm

**Example of labels for one code**
**REF 3810.15.010**

 D=14mm L=140mm  
 REVISION - Ti6Al4V


2022-08-31 QTY: 1 STERILE R Limacorporate S.p.a.

 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm  
**Uncemented Stem** EN  
**Stelo non Cementato** IT  
**Vastago no Cementado** ES  
**Zementlos Schaft** DE  
**Tige non Cimentée** FR  
**Haste não Cimentada** PT

**REF 3810.15.010**

07

**EN** Uncemented Stem  
**IT** Stelo non Cementato  
**ES** Vastago no Cementado  
**DE** Zementlos Schaft  
**FR** Tige non Cimentée  
**PT** Haste não Cimentada

D=14mm L=140mm 2022-08-31 Limacorporate S.p.A.

**REF 3810.15.010**

 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm

2022-08-31 LOT 2017 09554 STERILE R 1700266


 (20)07(240)93810150100(10)1709554(91)1700266(92)6826  
 (01)08033390003841(17)220831

07

**REF 3810.15.010**

 REVISION Hip/Anca - Ti6Al4V  
 Uncemented Stem/Stelo non Cementato  
 Dia=14mm L=140mm

2022-08-31 LOT 2017 09554 STERILE R 1700266


 (20)07(240)93810150100(10)1709554(91)1700266(92)6826  
 (01)08033390003841(17)220831

 Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy 0123


07

070

**WARNING**
**ENG:** (test according to ISO 7206-4)  
 Do not couple to lateralized necks

**ITA:** (test secondo ISO 7206-4)  
 Non accoppiabile a colli lateralizzanti

**DEU:** (Prüfung gemäß ISO7206-4)  
 Nicht geeignet für lateralisierte Halsteile

**FRA:** (test d'après ISO 7206-4)  
 Ne pas assembler aux cols latéralisants

**SPA:** (testado según ISO7206-4)  
 No acoplable a los cuellos lateralizados

**POR:** (testado segundo ISO7206-4)  
 No Adaptar a cósios lateralizados

**SECTION 2 – EXAMPLE OF LABELS FOR REVISION NECKS WITH LOCKING SCREW**

**Revision standard necks with locking screws**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




Code	Size
7515.15.005	H50
7515.15.010	H60
7515.15.020	H70
7515.15.030	H80
7515.15.040	H90
7515.15.050	H100
7515.15.060	H110

**Example of labels for one code**

REF 7515.15.010

12/14

REVISION - Ti6Al4V , UHMWPE



H60

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=60mm Taper 12/14

Neck + Screw **EN**

Collo + Vite **IT**

Cuello + Tornillo **ES**

Hals + Schraube **DE**

Col + Vis **FR**

Colo + Parafuso **PT**

2022-08-31 qty: 1 **STERILE EO** Limacorporate S.p.a.

REF 7515.15.010

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=60mm Taper 12/14

2022-08-31 **LOT** 2017 08292 **STERILE EO** 1700261



(20)07(240)97515150100(10)1708292(91)1700261(92)7473



(01)08033390017749(17)220831

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy **CE** 0123

REF 7515.15.010

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=60mm Taper 12/14

2022-08-31 **LOT** 2017 08292 **STERILE EO** 1700261



(20)07(240)97515150100(10)1708292(91)1700261(92)7473



(01)08033390017749(17)220831



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

**CE**  
0123



**Revision lateralized necks with locking screw**
**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).



**Code**

REF 7515.15.105

REVISION Hip - T16A4V, UHMWPE  
 Neck with Screw  
 H=50mm Lateralized-Taper 12/14

H50 Lat Size

2023-01-31 qty:1 STERILE EO Limacorporate S.p.a.

REF 7515.15.105

REVISION Hip - T16A4V, UHMWPE  
 Neck with Screw  
 H=50mm Lateralized-Taper 12/14

2023-01-31 LOT 2017 10114 STERILE EO 1800028

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy

04

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution Sterilized using Ethylene Oxide

IFU Consult instructions for use Manufacturer Sterilized using Irradiation

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy

CE 0123


Code	Size
7515.15.105	H50 Lat
7515.15.110	H60 Lat
7515.15.120	H70 Lat
7515.15.130	H80 Lat
7515.15.140	H90 Lat
7515.15.150	H100 Lat
7515.15.160	H110 Lat

**Example of labels for one code**

REF 7515.15.105

12/14

REVISION - Ti6Al4V - UHMWPE



H50  
Lat

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=50mm Lateralized-Taper 12/14

Neck + Screw **EN**

Collo + Vite **IT**

Cuello + Tornillo **ES**

Hals + Schraube **DE**

Col + Vis **FR**

Colo + Parafusc **PT**

2023-01-31 qty:1 **STERILEEO** Limacorporate S.p.a.

REF 7515.15.105

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=50mm Lateralized-Taper 12/14

2023-01-31 **LOT** 2017 10114 **STERILEEO** 1800028



(20)04(240)97515151050(10)1710114(91)1800028(92)4917



(01)08033390017800(17)230131

04

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy



REF 7515.15.105

REVISION Hip - Ti6Al4V , UHMWPE

Neck with Screw

H=50mm Lateralized-Taper 12/14

2023-01-31 **LOT** 2017 10114 **STERILEEO** 1800028



(20)04(240)97515151050(10)1710114(91)1800028(92)4917



(01)08033390017800(17)230131

04



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy




0123



**SECTION 3 – EXAMPLE OF LABELS FOR REVISION LR STEMS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




Code	Diameter	Length
3838.15.130	Dia 12	L 130
3840.15.130	Dia 14	L 130
3842.15.130	Dia 16	L 130
3838.15.165	Dia 12	L 165
3840.15.165	Dia 14	L 165
3842.15.165	Dia 16	L 165
3840.15.200	Dia 14	L 200
3842.15.200	Dia 16	L 200



**Example of label for one code**

REF 3838.15.130

REVISION LR-Ti6Al4V



Dia  
12  
L130

REVISION LR Hip/Anca Ti6Al4V  
Cemented Stem / Stelo Cementato  
Dia=12mm L=130mm  
Cemented Stem **EN**  
Stelo Cementato **IT**  
Vastago Cementado **ES**  
Zementschaft **DE**  
Tige Cimentée **FR**  
Haste Cimentada **PT**

2021-12-31 qty: 1 STERILE R Limacorporate S.p.a.

REF 3838.15.130

REVISION LR Hip/Anca Ti6Al4V  
Cemented Stem / Stelo Cementato  
Dia=12mm L=130mm

2021-12-31 LOT 2016 15822 STERILE R 1600344



(20)01(240)93838151300(10)1615822(9)1600344(92)6872  
(01)08033390045827(17)211231

REF 3838.15.130

REVISION LR Hip/Anca Ti6Al4V  
Cemented Stem / Stelo Cementato  
Dia=12mm L=130mm

2021-12-31 LOT 2016 15822 STERILE R 1600344



(20)01(240)93838151300(10)1615822(9)1600344(92)6872  
(01)08033390045827(17)211231

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123


Single use	Do not use if package is damaged	Keep away from sunlight	Temperature limit
Use By	Reference Number	Keep dry	
Lot Number	Caution	Sterilized using Ethylene Oxide	
IFU Consult instructions for use	Manufacturer	Sterilized using irradiation	

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

**SECTION 4 – EXAMPLE OF LABELS FOR REVISION LR MODULAR AUGMENTS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




Code	Size
3801.15.060	H60
3801.15.120	H120

Example of label for one code

**REF 3801.15.060**

REVISION LR-Ti6Al4V



**H60**

2021-10-31 **qty:** 1 **STERILE R** **Limacorporate S.p.a.**

REVISION LR Hip/Anca Ti6Al4V  
Modular Spacer/Distanziale Modulare  
H=60 - Additional/Aggiuntivo

**Modular Spacer** EN  
**Distanziale Modulare** IT  
**Distanciador Modular** ES  
**Modulunterlarge** DE  
**Entretoise Modulaire** FR  
**Distanciador Modular** PT

**REF 3801.15.060**

REVISION LR Hip/Anca Ti6Al4V  
Modular Spacer/Distanziale Modulare  
H=60 - Additional/Aggiuntivo

2021-10-31 **LOT** 2016 02905 **STERILE R** 1600293



**REF 3801.15.060**

REVISION LR Hip/Anca Ti6Al4V  
Modular Spacer/Distanziale Modulare  
H=60 - Additional/Aggiuntivo

2021-10-31 **LOT** 2016 02905 **STERILE R** 1600293



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy **CE** 0123

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult Instructions for use
- Manufacturer
- Sterilized using irradiation


Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy **CE** 0123



**SECTION 5 – EXAMPLE OF LABELS FOR MODULUS STEMS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




The colored band is added on the labels of Taper A stems and necks to identify the allowed combinations, as reported in the Surgical Technique.

Code	Diameter	Taper
4310.15.010	Dia 13	A
4310.15.020	Dia 14	A
4310.15.030	Dia 15	A
4310.15.045	Dia 16	B
4310.15.050	Dia 17	B
4310.15.060	Dia 18	B
4310.15.070	Dia 19	B
4310.15.080	Dia 20	B
4310.15.090	Dia 21	B
4310.15.100	Dia 22	B
4310.15.110	Dia 23	B
4310.15.120	Dia 24	B
4310.15.130	Dia 25	B
4310.15.140	Dia 26	B

Example of label for one code

REF 4310.15.010

MODULUS - Ti6Al4V



Dia.  
13  
A

MODULUS Hip System  
Stem / Stelo - Ti6Al4V

Dia=13mm - Taper / Cono A

Uncemented Stem **UN**

Stelo non Cementato **IT**

Vastago no Cementado **ES**

Zementlos Schaft **DE**

Tige non Cimentée **FR**

Haste não Cimentada **PT**

2022-03-31 qty: 1 STERILE R Limacorporate S.p.a.

REF 4310.15.010

MODULUS Hip System

Stem / Stelo - Ti6Al4V

Dia=13mm - Taper / Cono A

2022-03-31 LOT 2016 14390 STERILE R 1700095



(20)03(240)94310150100(10)1614390(91)1700095(92)5368

(01)08933390064749(17)228331



REF 4310.15.010

MODULUS Hip System

Stem / Stelo - Ti6Al4V

Dia=13mm - Taper / Cono A

2022-03-31 LOT 2016 14390 STERILE R 1700095



(20)03(240)94310150100(10)1614390(91)1700095(92)5368



(01)08933390064749(17)228331

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123



- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy


CE 0123



**SECTION 6 – EXAMPLE OF LABELS FOR MODULUS-R STEMS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




The colored band is added on the labels of stems with Ø up to 20mm to identify the allowed combinations, as reported in the Surgical Technique.

Code	Diameter	Height	Taper
4330.15.160	D16	H120	B
4330.15.170	D17	H120	B
4330.15.180	D18	H120	B
4330.15.190	D19	H140	B
4330.15.200	D20	H140	B
4330.15.210	D21	H140	B
4330.15.220	D22	H140	B
4330.15.230	D23	H140	B
4330.15.240	D24	H140	B
4330.15.250	D25	H140	B
4330.15.260	D26	H140	B

Example of labels for one code

REF 4330.15.160

MODULUS-R - Ti6Al4V



D16  
H120  
B

MODULUS-R Hip System  
Stem / Stelo - Ti6Al4V  
Dia=16mm H120mm – Taper / Cono B  
Uncemented Stem **EN**  
Stelo non Cementato **IT**  
Vastago no Cementado **ES**  
Zementlos Schaft **DE**  
Tige non Cimentée **FR**  
Haste não Cimentada **PT**

2022-03-31 qty: 1 STERILE R Limacorporate S.p.a.

REF 4330.15.160

MODULUS-R Hip System  
Stem / Stelo - Ti6Al4V  
Dia=16mm H120mm – Taper / Cono B

2022-03-31 LOT 2016 07417 STERILE R 1700080



(20)01(240)94330151600(10)1607417(91)1700080(92)5382



(01)08033390085311(17)220331

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123



REF 4330.15.160

MODULUS-R Hip System  
Stem / Stelo - Ti6Al4V

Dia=16mm H120mm – Taper / Cono B

2022-03-31 LOT 2016 07417 STERILE R 1700080



(20)01(240)94330151600(10)1607417(91)1700080(92)5382



(01)08033390085311(17)220331

01



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

0123



**SECTION 7 – EXAMPLE OF LABELS FOR MODULUS NECKS WITH LOCKING SCREW**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).




The colored band is added on the labels of Taper A stems and necks to identify the allowed combinations, as reported in the Surgical Technique.

Code	Size	CCD Angle	Taper
7590.15.010	S	125	A
7590.15.020	L	125	A
7590.15.030	S	125	B
7590.15.040	L	125	B
7590.15.050	S	125	B (Large neck)
7590.15.060	L	125	B (Large neck)
7595.15.010	S	135	A
7595.15.020	L	135	A
7595.15.030	S	135	B
7595.15.040	L	135	B
7595.15.050	S	135	B (Large neck)
7595.15.060	L	135	B (Large neck)



**Example of labels for one code**

**REF 7590.15.010**  
Taper 12/14  
MODULUS - Ti6Al4V, UHMWPE



**S**  
**125**  
**A**

2022-05-31 qty: 1 **STERILE** **EO** **Limacorporate S.p.a.**

MODULUS Hip System  
Neck / Collo - Ti6Al4V, UHMWPE  
S - 125 Deg./Gradi - Taper/Cono A

Neck + Screw **EN**  
Collo + Vite **IT**  
Cuello + Tornillo **ES**  
Hals + Schraube **DE**  
Col + Vis **FR**  
Colo + Parafuso **PT**

**REF 7590.15.010**   03

MODULUS Hip System  
Neck / Collo - Ti6Al4V, UHMWPE  
S - 125 Deg./Gradi - Taper/Cono A

2022-05-31 **LOT** 2016 16364 **STERILE** **EO** 1700171



201603240197590150100101616364(91)1700171(92)6939  
(01)08033390018197(17)220531

**Limacorporate S.p.A. - Via Nazionale, 52** **CE**  
33038 Villanova di San Daniele - Udine - Italy 0123

**REF 7590.15.010**  03  
MODULUS Hip System  
Neck / Collo - Ti6Al4V, UHMWPE  
S - 125 Deg./Gradi - Taper/Cono A

2022-05-31 **LOT** 2016 16364 **STERILE** **EO** 1700171



201603240197590150100101616364(91)1700171(92)6939  
(01)08033390018197(17)220531

 Single use	 Do not use if package is damaged	 Keep away from sunlight	 Temperature limit
 Use By	 Reference Number	 Keep dry	
 Lot Number	 Caution	 Sterilized using Ethylene Oxide	
 IFU Consult instructions for use	 Manufacturer	 Sterilized using irradiation	

**Limacorporate S.p.A. - Via Nazionale, 52** **CE**  
33038 Villanova di San Daniele - Udine - Italy 0123 