

2.2 Compatibles guide 0.018" (0,46 mm)

Sterling™
Ballon de dilatation

Monorail™

Références		Longueur Corps 135 cm	Longueur Corps 150 cm	Diamètre Ballon (mm)	Longueur Ballon (mm)	Compatibilité Introduceur (F)	Pression de Rupture (ATM) (kPa)	
Longueur Corps 80 cm	Longueur Corps 90 cm							
			H74939031202010*	2	20	4	14	1418
			H74939031204010*	2	40	4	14	1418
			H74939031206010*	2	60	4	14	1418
	H74939147208090*		H74939147208010*	2	80	4	14	1418
	H74939147201090*		H74939147201010*	2	100	4	14	1418
	H74939147201290*		H74939147201210*	2	120	4	14	1418
	H74939147201590*		H74939147201510*	2	150	4	14	1418
			H74939031202210*	2	220	4	14	1418
			H74939031252010*	2,5	20	4	14	1418
			H74939031254010*	2,5	40	4	14	1418
			H74939031256010*	2,5	60	4	14	1418
	H74939147258090*		H74939147258010*	2,5	80	4	14	1418
	H74939147251090*		H74939147251010*	2,5	100	4	14	1418
	H74939147251290*		H74939147251210*	2,5	120	4	14	1418
	H74939147251590*		H74939147251510*	2,5	150	4	14	1418
			H74939031252210*	2,5	220	4	14	1418
		H74939031302010		3	20	4	14	1418
		H74939031303010		3	30	4	14	1418
		H74939031304010		3	40	4	14	1418
		H74939031306010		3	60	4	14	1418
	H74939147308090*		H74939147308010*	3	80	4	14	1418
	H74939147301090*		H74939147301010*	3	100	4	14	1418
	H74939147301290*		H74939147301210*	3	120	4	14	1418
	H74939147301590*		H74939147301510*	3	150	4	14	1418
			H74939031302210*	3	220	4	14	1418
		H74939031352010		3,5	20	4	14	1418
		H74939031353010		3,5	30	4	14	1418
		H74939031354010		3,5	40	4	14	1418
			H74939031356010*	3,5	60	4	14	1418
	H74939147358090*		H74939147358010*	3,5	80	4	14	1418
	H74939147351090*		H74939147351010*	3,5	100	4	14	1418
	H74939147351290*		H74939147351210*	3,5	120	4	14	1418
	H74939147351590*		H74939147351510*	3,5	150	4	14	1418
			H74939031352210*	3,5	220	4	14	1418
H74939031401080				4	10	4	14	1418
H74939031401580				4	15	4	14	1418
H74939031402080		H74939031402010		4	20	4	14	1418
H74939031403080		H74939031403010		4	30	4	14	1418
H74939031404080		H74939031404010		4	40	4	14	1418
		H74939031406010		4	60	4	14	1418
	H74939147408090*		H74939147408010*	4	80	4	14	1418
	H74939147401090*		H74939147401010*	4	100	4	14	1418
	H74939147401290*		H74939147401210*	4	120	4	14	1418
	H74939147401590*		H74939147401510*	4	150	4	14	1418
			H74939031402210*	4	220	4	14	1418
H74939031451580				4,5	15	4	14	1418
H74939031452080		H74939031452010		4,5	20	4	14	1418
H74939031453080		H74939031453010		4,5	30	4	14	1418
H74939031454080		H74939031454010		4,5	40	4	14	1418
H74939031501080				5	10	4	14	1418
H74939031501580				5	15	4	14	1418
H74939031502080		H74939031502010		5	20	4	14	1418
H74939031503080		H74939031503010		5	30	4	14	1418
H74939031504080		H74939031504010		5	40	4	14	1418
		H74939031506010		5	60	4	14	1418

*Sterling™ SL

Sterling™ Suite

Ballon de dilatation

Monorail™

Longueur Corps 80 cm	Références			Diamètre Ballon (mm)	Longueur Ballon (mm)	Compatibilité Introducteur (F)	Pression de Rupture (ATM) (kPa)	
	Longueur Corps 90 cm	Longueur Corps 135 cm	Longueur Corps 150 cm					
H74939031551580				5,5	15	4	14	1418
H74939031552080		H74939031552010		5,5	20	4	14	1418
H74939031553080		H74939031553010		5,5	30	4	14	1418
H74939031554080		H74939031554010		5,5	40	4	14	1418
H74939031601080				6	10	4	14	1418
H74939031601580				6	15	4	14	1418
H74939031602080		H74939031602010		6	20	4	14	1418
H74939031603080		H74939031603010		6	30	4	14	1418
H74939031604080		H74939031604010		6	40	4	14	1418
		H74939031606010		6	60	4	14	1418
H74939031651580				6,5	15	4	14	1418
H74939031652080		H74939031652010		6,5	20	4	14	1418
H74939031653080		H74939031653010		6,5	30	4	14	1418
H74939031654080		H74939031654010		6,5	40	4	14	1418
H74939031701080				7	10	4	14	1418
H74939031701580				7	15	4	14	1418
H74939031702080		H74939031702010		7	20	4	14	1418
H74939031703080		H74939031703010		7	30	4	14	1418
H74939031704080		H74939031704010		7	40	4	14	1418
		H74939031706010		7	60	5	14	1418
H74939031801080				8	10	4	12	1216
H74939031801580				8	15	4	12	1216
H74939031802080		H74939031802010		8	20	4	12	1216
H74939031803080		H74939031803010		8	30	4	12	1216
H74939031804080		H74939031804010		8	40	4	12	1216
		H74939031806010		8	60	5	12	1216