

**pITA<sup>®</sup>**

Rapid Exchange  
Neurovascular PTA Balloon Catheter

**PITA-125-10**

Rapid Exchange Neurovascular PTA Balloon Catheter  
Neurovaskulärer Rapid Exchange Ballonkatheter  
Cathéter à ballonnet pour ATP neurovasculaire à échange rapide  
Catetere da PTA a palloncino Rapid Exchange neurovascolare  
Catéter con balón neurovascular de intercambio rápido para ATP  
Cateter Balão para PTA Neurovascular de Troca Rápida  
нейроваскулярный баллонный катетер быстрой замены для ЧТА

**phenox GmbH**

Lise-Meitner-Allee 31 - 44801 Bochum, Germany  
Tel: +49 (0) 234 / 36919 - 0 Fax: +49 (0) 234 36919 - 19  
www.phenox.net

**1.25 mm**


Balloon Diameter

**10 mm**

Balloon Length

**160 cm**

Catheter Length



| A:<br>Balloon Diameter<br>[mm] | B:<br>Balloon Length<br>[mm] | C:<br>Catheter Length<br>[mm] | D:<br>Distal Shaft<br>Diameter | E:<br>Proximal Shaft<br>Diameter |
|--------------------------------|------------------------------|-------------------------------|--------------------------------|----------------------------------|
| 1.25                           | 10                           | 1600                          | 2.5F / 0.83 mm                 | 1.9F / 0.63 mm                   |

|                                                |                      |
|------------------------------------------------|----------------------|
| Guidewire                                      | 0.014"<br>/ 0.36 mm  |
| Guiding Catheter<br>Compatibility<br>(Min. ID) | 0.055"<br>/ 1.4 mm   |
| Rated Burst<br>Pressure (RBP)                  | 18 atm<br>/ 1824 kPa |
| Nominal Pressure<br>(NOM)                      | 10 atm<br>/ 1013 kPa |

|     |      |      |      |      |      |      | NOM  |      |      |      |      |      |      |      |      |      |      | RBP |  |  |
|-----|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|-----|--|--|
| atm | 4    | 5    | 6    | 7    | 8    | 9    | 10   | 11   | 12   | 13   | 14   | 15   | 16   | 17   | 18   | 19   | 20   |     |  |  |
| mm  | 1.15 | 1.17 | 1.19 | 1.21 | 1.22 | 1.24 | 1.25 | 1.26 | 1.27 | 1.27 | 1.28 | 1.29 | 1.29 | 1.30 | 1.30 | 1.31 | 1.32 |     |  |  |

**REF** PITA-125-10

**LOT** 15118809

**2017-09**

**Contents**  
1x pITA  
2x Clip  
1x Cannula

pITA **REF** PITA-125-10 **LOT** 15118809



(01)04260123781118(17)170918(10)15118809

pITA **REF** PITA-125-10 **LOT** 15118809



(01)04260123781118(17)170918(10)15118809

pITA **REF** PITA-125-10 **LOT** 15118809



(01)04260123781118(17)170918(10)15118809

LAB-B837-PITA Rev.D

**pITA<sup>®</sup>** Rapid Exchange  
Neurovascular PTA Balloon Catheter

1.25mm x 10 mm



**REF** PITA-125-10

**LOT** 15118809 **2017-09**

**CE0297**

**STERILE**  
Sterilized with Ethylene  
Oxide Gas

**Consult**  
instruction for use

**Caution**

**Do not reuse**

**Keep dry**

**Keep away**  
from sunlight

**Do not use if the**  
packaging is damaged

**Do not**  
resterilize

**Non-pyrogenic**

**+10°C**

**+40°C**

