

FoxCross

PTA Catheter



(01)0000000000000(17)010131(10)000000(91)3000

REF LV_ITNR

LOT *VOID

LV_EXPD

BALLOON

3.0 mm

BALLOON

20 mm

USAB CAT

80 cm

SHEATH

5 Fr
1.67 mm

FoxCross


PTA Catheter

CONTENTS 1 FoxCross OTW PTA Catheter

REF LV_ITNR

LOT *VOID

LV_EXPD



BALLOON

3.0 mm

BALLOON

20 mm

RBP

18 bar
1800 kPa

USAB CAT

80 cm

SHEATH

5 Fr
1.67 mm

REC GW

0.035"
0.89 mm

STERILE R

NON PYROGENIC



15°C 30°C



DIST

Abbott Vascular Devices Holland B.V.,
Maastricht Airport, Bamfordweg 1,
NL-6235 NS Uiterstraten, The Netherlands
Phone: +31 43 358 6751
FAX: +31 43 358 6759Abbott Laboratories Vascular Enterprises Ltd.,
Dublin, Beringen Branch, Anthoptstrasse,
CH-8222 Beringen, Switzerland,
Phone: +41 52 674 8550
FAX: +41 52 674 8540**Abbott**

© 2008 Abbott Laboratories

LA 010920001 001

DO NOT EXCEED!

NPR RBP **

Balloon Compliance ϕ 3.0			
	Pressure		Diameter [mm]
	[bar]	[kPa]	
	4	400	2.6
	5	500	2.7
	6	600	2.8
	7	700	2.8
	8	800	2.9
	9	900	3.0
*	10	1000	3.0
	11	1100	3.0
	12	1200	3.1
	13	1300	3.1
	14	1400	3.2
	15	1500	3.2
	16	1600	3.2
	17	1700	3.3
**	18	1800	3.3
	19	1900	3.3
	20	2000	3.4
	21	2100	3.4
	22	2200	3.4
	23	2300	-
	24	2400	-

FoxCross

PTA Catheter

 ϕ 3.0 mm

k 20 mm

FoxCross

PTA Catheter

 ϕ 3.0 mm

k 20 mm

FoxCross

PTA Catheter

BALLOON 3.0 mm

BALLOON 20 mm

ALIVE Ltd, Dublin, Beringen Branch, CH-8222 Beringen, Phone: +41 52 674 8550 FAX: +41 52 674 8540



(01)00000000000000



(17)010131(10)000000(91)3000



(01)00000000000000



(17)010131(10)000000(91)3000