


PRODUCT CODE	SC60
Description	Jeweller Style Forceps

Component	Material of the device	Stainless steel
	Material in contact with the patient	Stainless steel
Destination	Dimensions	0.12 mm fine pointed ends. 32.26mm straight jaws. Handle length: 82mm. Overall Length: 115mm
	Intended use	Used to hold ocular tissue, remove eye lashes, for suture tying and to remove a 6mm window in the anterior chamber and to remove part of an artificial intraocular lens
Sterilization	Single use or reusable?	Single use
	Sterile or not sterile?	Sterile
	If sterile : sterilization procedure	Gamma irradiated
	If sterile : reusable?	No
	If non sterile : process of sterilization	N/A
Packaging and storage	Primary packaging component	Peel pouch (Tyvek / Pet-PE film)
	Quantity per sterile unit	1
	Clinical box component	cardboard box
	Quantity per clinical box	10
	Storage conditions	Humidity 20-80% Temp -5°C/50°C
	Shelf life	3 Years
Regulatory	Legal Manufacturer	SURGITRAC SIHONG Ltd
	EC certificate	EC Certificate Directive 93/42/EEC Annex V : DD60147761 0001
	Class of the device (according 93/42/EEC MDD)	Class IIa
	Classification Rule (according 93/42/EEC MDD)	6
	GMDN Code	62489
	Presence of biological tissues or substances of animal origin?	NO biological tissues or substances of animal origin
	Latex or Latex free?	Latex free
	Presence of pharmaceuticals/Medicines	No pharmaceutical or medicinal substance



Date : 12/5/2025

Validated by : Marie SUET Quality and Regulatory Affairs Manager
quality@surgitrac-europe.com
 2 rue Hélène Boucher – 35235 Thorigné-Fouillard, France.